



NCTTRAC Distribution & Implementation Plan – Pediatric Readiness Training Equipment

I. Purpose

- A. Provide pediatric Quality Cardiopulmonary Resuscitation (QCPR) training manikins to support the Pediatric Readiness Project.
- B. Establish procedures for the distribution, management, and disposition of the Pediatric Readiness Training Equipment purchased and distributed by NCTTRAC.

II. Description / Distribution of Asset

- A. Distribution of the training equipment will be facilitated at the NCTTRAC office. The Point of Contact (POC) information will be available on the NCTTRAC website.
- B. NCTTRAC member organizations may reserve (subject to availability) up to **two** pediatric QCPR manikins for training at their organization for up to **two** weeks.
- C. Members will be responsible for:
 - Providing NCTTRAC with the required contact information
 - Maintaining equipment integrity and cleanliness
 - Facilitating Pediatric Readiness Simulation training
 - Returning equipment to the NCTTRAC office within the **two-week** requirement
 - Notifying NCTTRAC staff of damaged or malfunctioning equipment within 48 hours of identifying the issue in writing
- D. Asset management will be facilitated by utilizing an electronic check-in/out process.
- E. The *Equipment Damages Liability Agreement* must be signed by the requesting organization prior to equipment checkout.
- F. Each Resusci Junior QCPR manikin includes one gray carrying bag and the following equipment:
 - Junior Resusci QCPR manikin with jacket and shorts(1)
 - Pulse bulb (1)
 - Important Product Information (1)
 - User Guide (1)
 - QCPR App Instructions (1)
 - USB-C power charging cord (1)
 - AA size batteries (4)
- G. Each Resusci Baby QCPR manikin includes one gray plastic carrying case and the following equipment:
 - Full-Body manikin dressed in a onesie
 - QCPR App Instructions (1)
 - Pulse bulb (1)
 - Disposable Lungs (2)
 - Resusci Baby Blanket (1)
 - QCPR armband (1)
 - Important Product Information (1)
 - User Guide (1)
 - Airway Assemblies (2)
 - USB-C power charging cord (1)

III. Asset Disposition

NCTTRAC staff will consult committee members regarding the replacement or demobilization of a manikin should it become undeployable.

IV. Limitations to Possession of Asset

The pediatric QCPR manikins are intended for use by healthcare providers in TSA-E within the scope and limits of this plan. Neither NCTTRAC nor this plan certifies or authorizes an individual in Pediatric Readiness. Individuals or organizations who seek Pediatric Readiness training can request more information through the NCTTRAC ticket system at www.ncttrac.org.